

# MINDACT Trial Results

(Microarray In Node negative and 1-3 positive lymph node Disease may Avoid Chemotherapy)


MINDACT is the first and only Phase III prospective, randomized controlled clinical trial comparing a genomic breast cancer recurrence assay to clinicopathological risk assessment (standard of care).

## Baseline characteristics of the MINDACT study population

	LOW RISK (no chemo)* C=↓ MP=↓ n = 2,745 (41%)	OVERTREATMENT GROUP C=↑ MP=↓ CPFs HIGH MammaPrint LOW n = 1,550 (23%)	UNDERTREATMENT GROUP C=↓ MP=↑ CPFs LOW MammaPrint HIGH n = 592 (9%)	HIGH RISK (chemo)† C=↑ MP=↑ n = 1,806 (27%)
IHC Subtype	96% Luminal • HER2+ : 4%	91% Luminal • HER2+ : 8% • Triple Neg : 1%	79% Luminal • HER2+ : 12% • Triple Neg : 9%	50% Luminal • HER2+ : 19% • Triple Neg : 31%
Grade	Grade 1 or 2 98%	Grade 2 or 3 94%	Grade 1 or 2 85%	Grade 3 76%
Node Positive	6%	48%	2%	26%
Tumor Size >2cm	4%	58%	2%	48%
5-year DMFS Distant Metastasis Free Survival	97.6% (96.9 - 98.1) - no ACT	94.7% (92.5 - 96.2) - no ACT	93.9% (89.6 - 96.5) - no ACT	90.6% (89.0 - 92.0) - ACT

Primary Analysis Population

## Data from the Primary Analysis: 5-year DMFS for patients with high risk clinical features identified as Low Risk by MammaPrint, without chemotherapy




# MINDACT Trial Results

(Microarray In Node negative and 1-3 positive lymph node Disease may Avoid Chemotherapy)

## Efficacy of chemotherapy was determined by comparing DMFS at 5 years (All patients cHigh / MP Low)

**OVERTREATMENT GROUP**

C=↑ MP=↓  
CPFs HIGH  
MammaPrint LOW



- MINDACT data shows that in the presence of high risk clinicopathologic features, there is **no significant or clinically meaningful benefit to adding chemotherapy for MammaPrint Low Risk patients.**

## MINDACT Results Summary

MammaPrint was developed independent of clinicopathologic factors so it's the only assay that provides added value to the treatment decision and can be used in conjunction with clinical factors.

MINDACT demonstrates that MammaPrint is the only breast cancer genomic assay that selectively identifies early-stage breast cancer patients with high risk features (Lymph node positivity (LN 1-3), grade II or III, size >2cm) as Low Risk with no clinically meaningful benefit of chemotherapy.

MINDACT proves that MammaPrint is superior as a standalone test to the current standard of care in predicting the benefit of chemotherapy for early-stage breast cancer patients, minimizing the risk of undertreatment.



MINDACT provides indisputable Level 1A, phase III, randomized prospective clinical evidence that proves the superiority of MammaPrint 70-Gene Assay in predicting the benefit of chemotherapy in early-stage breast cancer patients when compared to clinicopathological risk assessment.

### References

- Piccart M. *AACR Podium Presentation*, April 18th, 2016